



# Recombinant Escherichia coli O17:K52:H18 Glycogen debranching enzyme (glgX)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP001446EOG   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | B7NE41  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | <p>MTQLAIGKPA PLGAHYDGQG VNFTLFSAHA ERVELCIFDA NGQEHRYDLP<br/> GHSGLDIWHGY LPDARPGRLY GYRVHGPWQP AEGHRFNPAK LLIDPCARQI<br/> DGEFKDNPLL HAGHNEPDYR DNAAIAPKCV VVDHYDWD DAPRTPWGS<br/> TIIYEAHVKG LTYLHPEIPV EIRGTYKALG HPVMINYLKQ LGITALELLP<br/> VAQFASEPRL QRMGLSNIYWG YNPVAMFALH PAYACSPETA LDEFRDAIKA<br/> LHKAGIEVIL DIVLNHSAEL DLDGPLFSLR GIDNRSYYWI REDGDYHNWT<br/> GCGNTLNLSH PAVVDYASAC LRYWVETCHV DGFRFDLAAV MGRTPEFRQD<br/> APLFTAIQNC PVLSQVKLIA EPWDIAPGGY QVGNFPPLFA EWNDHFRDAA<br/> RRFWLHYDLP LGAFAGRFAA SSDVFKRNGR LPSAAILNLT AHDGFTLRDC<br/> VCFNHKHNEA NGEENRDGTN NNYSNNHGKE GLGGTLDLVE RRRDSIHALL<br/> TTLLSQGTP MLLAGDEHGH SQHGNNNAYC QDNQLTWLDW<br/> SQASSGLTAF TAALHLRKR IPALVENRWW EEGDGNVRWL NRYAQPLSTD<br/> EWQNGPKQLQ ILLSDRFLIA INATLEVTEI VLPAGEWHAI PPFAGEDNPV<br/> ITAVWQGPAH GLCVFQR</p> |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | glgX  |
| <b>Protein Names</b>     | Recommended name: Glycogen debranching enzyme EC= 3.2.1.- Alternative name(s): Glycogen operon protein GlgX   |
| <b>Expression Region</b> | 1-657   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.