



Recombinant Escherichia coli O17:K52:H18 N-acetyl-D-glucosamine kinase (nagK)

Product Code	CSB-MP015415EOG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7NAZ7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MYYGFDIGGT KIALGVFDSG RQLQWEKRVP TPRDSYDAFL DAVCELVAEA DQRFGCKGSV GIGIPGMPET EDGTLYAANV PAASGKPLRA DLSARLDRDV RLDNDANCFALSEAWDDEFT QYPLVMGLIL GTGVGGGLIF NGKPITGKSY ITGEFGHMRL PVDALTMMLDFPLRRCGCG QHGCNIENYLS GRGFAWLYQH YYHQPLQAPE IIALYDQGDE QARAHVERYL DLLAVCLGNI LTIVDPDLVV IGGGLSNFPA ITTQLADRLP RHLLPVARVP RIERARHGDA GGMRGAAFLH LTD
Source	Mammalian cell
Target Names	nagK
Protein Names	Recommended name: N-acetyl-D-glucosamine kinase EC= 2.7.1.59 Alternative name(s): GlcNAc kinase
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.