



# Recombinant Escherichia coli O127:H6 Glycogen debranching enzyme (glgX)

<b>Product Code</b>	CSB-EP001446EOB-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7UKY8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MTQLAIGKPT PLGAHYDGQG VNFTLFSAHA ERVELCVFDA NGQEHRYDLP  GHSGLDIWHGY LPDARPGRLY GYRVHGPWQP AEGHRFNPAK LLIDPCARQI  DGEFKDNPLL HAGHNEPDYR DNAAIAPKCV VVDHYDWD DAPRTPWGS  TIIYEAHVKG LTYLHPEIPV EIRGTYKALG HPVMINYLKQ LGITALELLP  VAQFASEPRL QRMGLSNIYWG YNPVAMFALH PAYACSPETA LDEFRDAIKA  LHKAGIEVIL DIVLNHSAEL DLDGPLFSLR GIDNRSYYWI REDGDYHNWT  GCGNTLNLSH PAVVDYASAC LRYWVETCHV DGFRFDLAAV MGRTPEFRQD  APLFTAIQNC PVLSQVKLIA EPWDIAPGGY QVGNFPPLFA EWNDHFRDAA  RRFWLHYDLP LGVFAGRFAA SSDVFKRNGR LPSAAILNLT AHDGFTLRDC  VCFNHKHNEA NGEENRDGTN NNYSNNHGKE GLGGTLDLVE RRRDSIHALL  TTLSSQGTP MLLAGDEHGH SQHGNNNAYC QDNQLTWLDW  SQASSGLTAF TAALIHRLKR IPALMENRWW EEGDGNVRWL NRYAQLSTD  EWQNGPKQLQ ILLSDRFLIA INATLEVTEI VLPAGEWHAI PPFAGEDNPV  ITAVWQGPAH GLCVFQR</p>
<b>Source</b>	E.coli
<b>Target Names</b>	glgX
<b>Protein Names</b>	Recommended name: Glycogen debranching enzyme EC= 3.2.1.- Alternative name(s): Glycogen operon protein GlgX
<b>Expression Region</b>	1-657
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.