



Recombinant Escherichia fergusonii Argininosuccinate synthase (argG)

Product Code	CSB-YP002234EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LR40
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MTTILKHLPV GQRIGIAFSG GLDTSAALLW MRQKGAVPYA YTANLGQPDE EDYDAIPRRA MEYGAENARL IDCRKQLVAE GIAAIQCGAF HNTTGGLTYF NTTPLGRAVT GTMLVAAMKE DGVNIWGDGS TYKGN DIERF YRYGLLTNAE LQIYKPWLDT DFIDELGGRH EMSEFMTACG FDYKMSVEKA YSTDSNMLGA THEAKDLEYL NSSVKIVNPI MGVKFWDES V KIPAE ETVR FEQGHPVALN GKTFSDDVEM MLEANRIGGR HGLGMSDQIE NRIIEAKSRG IYEAPGMALL HIAYERLLTG IHNEDTIEQY HAHGRQLGRL LYQGRWFDSQ ALMLRDSLQR WVASQITGEV TLELRGNDY SILNTISNNL TYKPERLTME KGDSVFSPDD RIGQLTMRNL DITDTREKLF GYAKTGLLSS SAASGVPQVE NLENKAK
Source	Yeast
Target Names	argG
Protein Names	Recommended name: Argininosuccinate synthase EC= 6.3.4.5 Alternative name(s): Citrulline--aspartate ligase
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.