



# Recombinant Escherichia coli O7:K1 tRNA dimethylallyltransferase (miaA)

<b>Product Code</b>	CSB-EP024529EON
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7NTM7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDISKASLP KAIFLMGPTA SGKTALAIEL RKILPVELIS VDSALIYKGM DIGTAKPNAE ELLAAPHRLD DIRDPSQAYS AADFRRDALA EMADITAAGR IPLLVGGTML YFKALLEGLS PLPSADPEVR ARIEQQAEEQ GWESLHRQLQ EIDPVAARI HPNDPQRLSR ALEVFFISGK TLTELTQTSG DALPYQVHQF AIAPASRELL HQRIEQRFHQ MLASGFEEV RALFARGDLH TDLPSIRCVG YRQMWSYLEG EISYDEMYYR GVCATRQLAK RQITWLRGWE GVHWLDSEKP EQARDEVLQV VGAIAG
<b>Source</b>	E.coli
<b>Target Names</b>	miaA
<b>Protein Names</b>	Recommended name: tRNA dimethylallyltransferase EC= 2.5.1.75 Alternative name(s): Dimethylallyl diphosphate:tRNA dimethylallyltransferase Short name= DMAPP:tRNA dimethylallyltransferase Short name= DMATase Isopentenyl-diph
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.