



Recombinant Escherichia fergusonii Porphobilinogen deaminase (hemC)

Product Code	CSB-MP010524EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LU53
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MLDNVLR IAT RQSPLALWQA HYVKDALMAK HPGLTVELVP MVTRGDVILD TPLAKVGGKG LFKVELEVAL LENRADI AVH SMKDVPEFEP QGLGLVTICE REDPRDAFVS NKYASLDEL P AGSIVGTSSL RRQCQLAERR PDLIIRSLRG NVGTRLSKLD NGEYDAIILA VAGLKRLGLE SRIRDALPPE VSLPAVGQGA VGIECRLDDT RTRELLAALN HPETALRVTA ERAMNTRLEG GCQVPIGSYA ELIGGEIWL R ALVGAPDGSQ MIRGERRGSP QDAEKMGISL AEELLNNGAR AILADVNGD APV
Source	Mammalian cell
Target Names	hemC
Protein Names	Recommended name: Porphobilinogen deaminase Short name= PBG EC= 2.5.1.61 Alternative name(s): Hydroxymethylbilane synthase Short name= HMBS Pre-uroporphyrinogen synthase
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.