



Recombinant Escherichia fergusonii Triosephosphate isomerase (tpiA)

Product Code	CSB-BP024102EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LVC6
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MRHPLVMGNW KLNGSRHMVH ELVSNLRKEL AGVAGCAVAI APPEMYIDMA KREAEGSHIM LGAQNVDLNL SGAFTGETSA AMLKDIGAQY IIIGHSERRT YHKESDELIA KFAVLKEQG LTPVLCIGET EAENEAGKTE EVCARQIDAV LKTQGAAAFE GAVIAYEPVW AIGTGKSATP AQAQAVHKFI RDHIAKVDAN IAEQVIIQYG GSVNASNAAE LFAQPDIDGA LVGGASLKAD AFAVIVKAAE AAKQA
Source	Baculovirus
Target Names	tpiA
Protein Names	Recommended name: Triosephosphate isomerase Short name= TIM EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.