



Recombinant Escherichia coli O127:H6 NAD-dependent malic enzyme (maeA)

Product Code	CSB-BP013633EOB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7URL9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MEPKTKKQRS LYIPYAGPVL LEFPLLKGS AFSMEERRNF NLLGLLPEVV ETIEEQAERA WIQYQGFKTE IDKHIYLRNI QDTNETLFYR LVNNHLDEMM PVIYTPTVGA ACERFSEIYR RSRGVFISYQ NRHNMDLILQ NVPNHNIKVI VVTDGERILG LGDQGIGGMG IPIGKLSLYT ACGGISPAYT LPVVLVDVGTN NQQLLNDPLY MGWRNPRITD DEYYEFVDEF IQAVKQRWPD VLLQFEDFAQ KNAMPLLNRY RNEICSFNDD IQGTAAVTVG TLIAASRAAG GQLSEKKIVF LGAGSAGCGI AEMIIAQTR EGLSEEAARQ KVFMVDRFGL LTDKMPNLLP FQTKLVQKRE NLSDWDTSD VLSLLDVVRN VKPDILIGVS GQTGLFTEEI IREMHKHCPR PIVMPLSNPT SRVEATPQDI IAWTEGNALV ATGSPFNPVV WKDKIYPIAQ CNAFIFPGI GLGVIASGAS RITDEMLMSA SETLAQYSPL VLNGEGLVLP ELKDIQKVS R AIAFAVGKMA QQQGVAVKTS AEALQQAIDD NFWHAEYRDY RRTSI
Source	Baculovirus
Target Names	maeA
Protein Names	Recommended name: NAD-dependent malic enzyme Short name= NAD-ME EC= 1.1.1.38
Expression Region	1-565
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.