



Recombinant Escherichia fergusonii Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LLU4
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MLNSFKLSLQ YILPKLWLTR LAGWGASKRA GWLTKLVIDL FVKYYKVDMMK EAQKPDTASY RTFNEFFVRP LRDEVRPIDT DPNVLVMPAD GVISQLGKIE EDKILQAKGH NYSLEALLAG NYLMADLFRN GTFVTTYLSP RDYHRVHMP NGILREMIYV PGDLFSVNHL TAQNVNLFARNERVICLFD TFGPMAQIL VGATIVGSIE TVWAGTITPP REGIIKRWTW PAGENDGSVA LLKGQEMGRF KLG
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.