



Recombinant *Desulfurococcus kamchatkensis* Cell division control protein 6 homolog (cdc6)

Product Code	CSB-BP005022DJT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B8D6H2
Product Type	Recombinant Protein
Immunogen Species	<i>Desulfurococcus amylolyticus</i> (strain DSM 18924 / JCM 16383 / VKM B-2413 / 1221n) (<i>Desulfurococcus kamchatkensis</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MADIIDELIR KRGFASRIFR NREVLHPDYI PDTLPHRENE IRRLAHELLV SAQGMRPSNV LIYGLTGTGK TVVVKYVVS LKEKASSLNK RLDYAYVNTR KLDTTYRVIA SIAQSIGLRV PHTGLAISEV YRRYINALDS WGGLHIIVLD EVDYYVKREG DDLIYKLVRA NEELSKAKIV LIGITNDVNF VENLDPRVRS SMGEIEMVFP PYNAEQLFTI LKQRAELAFN QGVIEDGVIS YCSALAAREH GDARRALDLL RVSGEIAERE GAERTIEYV KKATLEIEEG RIYQSVVTLTP LHQKLVLKKI VELVEAKEST TTGEVYTAYS NIMRNLYKYP LSLRRISEII SQLDMMGLII SEVVNRGKYG ITRIIKIRKD MLPVIKDALK DVEA
Source	Baculovirus
Target Names	cdc6
Protein Names	Recommended name: Cell division control protein 6 homolog Short name= CDC6 homolog
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.