



Recombinant Rat Kinesin heavy chain isoform 5C (Kif5c)

Product Code	CSB-MP012343RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56536
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MADPAECSIK VMCRFRPLNE AEILRGDKFI PKFKGEETVV IGQGKPYVFD RVLPPNTTQE QVYNACAKQI VKDVLEGYNG TIFAYGQTSS GKHTMEGKL HDPQLMGIIP RIAHDIFDHI YSM DENLEFH IKVSYFEIYL DKIRDLLDVS KTNLAVHEDK NRVPYVKGCT ERFVSSPEEV MDVIDEGKAN RHVAVTNMNE HSSRSHSIFL INIKQENVET EK KLSGKLYL VDLAGSEKV
Source	Mammalian cell
Target Names	Kif5c
Protein Names	Recommended name: Kinesin heavy chain isoform 5C Alternative name(s): Kinesin heavy chain neuron-specific 2
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.