



Recombinant Mouse CAP-Gly domain-containing linker protein 3 (Clip3)

Product Code	CSB-YP005555MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B9EHT4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MTKTDPAAPMA PPPRGEEEEEE EEEDEPVPEA PSPTQERRQK PVVHPSAPAP LPKDYAFTFF DPNDPACQEI LFDPKTTIPE LFAIVRQWVP QVQHKIDVIG NEILRRGCHV NDRDGLTDMT LLHYACKAGA HGVGDPAAAV RLSQQLLAG ADVTLSRWT NMNALHYAAY FDVPDLVRVL LKGARPRVVN STCSDFNHGS ALHIAASNLC LGAAKCLLEH GANPALRNRK GQVPAEVVPD PMDMSLDKAE AALVAKELRT LLEEAVPLSC TLPKVTLPNY DNVPGNLMLS ALGLRLGDRV LLDGQKTGTL RFCGTTEFAS GQWVGVELDE PEGKNDGSVG GVRFYICPPK QGLFASVSKV SKAVDAPPSS VTSTPRTPRM DFSRVTGKGR REHKGKKKSP SSPSLGLSQQ REGAKAEVGD QVLVAGQKQG IFRFYGKTDF APGYWYGIEL DQPTGKHDGS VFGVRYFTCA PRHGVFAPAS RIQRIGGSTD PPGDSVGAKK VHQVTMTQPK RTFTTVRTPK DIASENSISR LLFCCWFPWM LRAEMQS
Source	Yeast
Target Names	Clip3
Protein Names	Recommended name: CAP-Gly domain-containing linker protein 3 Alternative name(s): Cytoplasmic linker protein 170-related 59 kDa protein Short name= CLIP-170-related 59 kDa protein Short name= CLIPR-59
Expression Region	1-547
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.