



# Recombinant Human Poly (rC)-binding protein 4 (PCBP4)

<b>Product Code</b>	CSB-MP017521HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P57723
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSGSDGGLEE EPELSITLTL RMLMHGKEVG SIIGKKGETV KRIREQSSAR ITISEGSCPE RITTITGSTA AVFHAVSMIA FKLDEDLCAA PANGGNVSRP PVTLRLVIPA SQCGSLIGKA GTKIKEIRET TGAQVQVAGD LLPNSTERAV TVSGVPDAII LCVRQICAVI LESPPKGATI PYHPSLSLGT VLLSANQGF VQGQYGAVTP AEVTKLQQLS SHAVPFATPS VVPGLDPGTQ TSSQEFLVPN DLIGCVIGRQ GSKISEIRQM SGAHIKIGNQ AEGAGERHVT ITGSPVSIAL AQYLITACLE TAKSTSGGTP SSAPADLPAP FSPPLTALPT APPGLLGTPY AISLSNFIGL KPMPFLALPP ASPGPPPGLA AYTAKMAAAN GSKKAERQKF SPY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PCBP4
<b>Protein Names</b>	Recommended name: Poly(rC)-binding protein 4 Alternative name(s): Alpha-CP4
<b>Expression Region</b>	1-403
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a member of the KH-domain protein subfamily. Proteins of this subfamily, also referred to as alpha-CPs, bind to RNA with a specificity for C-rich pyrimidine regions. Alpha-CPs play important roles in post-transcriptional activities and have different cellular distributions. This gene is induced by the p53 tumor suppressor, and the encoded protein can suppress cell proliferation by inducing apoptosis and cell cycle arrest in G(2)-M. This gene's protein is found in the cytoplasm, yet it lacks the nuclear localization signals found in other subfamily members. Multiple alternatively spliced transcript variants have been described, but the full-length nature for only some has been determined.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.