



Recombinant Pseudomonas fluorescens Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033FFU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C3KDM6
Product Type	Recombinant Protein
Immunogen Species	Pseudomonas fluorescens (strain SBW25)
Purity	>85% (SDS-PAGE)
Sequence	MKKQLFILSQ YLLPHHLLSR LAGCIAECRV RWFKNAFTAW FAKRYQVDMS QALVEDLTAY EHFNAFFTRA LKDGARPLDQ TPGAVLSPAD GAVSQLGPIE HGRVFAQKGH SFSVLELLGG DAALSAPFMG GDFATVYLSP KDYHRVHMPL AGTLREMYI PGRIFSVNQT TAENVPELFA RNERVACIFD TERGPMVVL VGAMIVASIE TVWAGLVTPP KRELKTFRYD EAARAPIHLE KGAE LGRFKL G
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.