



# Recombinant *Hamiltonella defensa* subsp. *Acyrtosiphon pisum* Phosphatidylserine decarboxylase proenzyme (psd)

<b>Product Code</b>	CSB-MP018033HUG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	C4K3T9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Hamiltonella defensa</i> subsp. <i>Acyrtosiphon pisum</i> (strain 5AT)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLKCLKIKLQ YLLPQHALTR LAGWAANKKA GWLTQQVIKI FASYRVNIK EAQYTEFSAY SSFNEFFIRS LRSDVRPIAA GDNILVQPAD GTISESGRID KDKILQAKGH DYSLEALLAG QSLLAKEFEN GQFVTTYLAP GDYHRVHMPC DGVLREMIYV PGSLFSVNAF FAENVPYLFARNERVICIFN TAFGTMAQIL VGAMIVGSIE TLWSGPVDSQ RKGIIQRWVY PDEGSEAPII LKKGEEEMGLF KLG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	psd
<b>Protein Names</b>	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
<b>Expression Region</b>	1-253
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.