



Recombinant Tolumonas auensis Lipoyl synthase (lipA)

Product Code	CSB-YP012927TOQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C4LCC5
Product Type	Recombinant Protein
Immunogen Species	Tolumonas auensis (strain DSM 9187 / TA4)
Purity	>85% (SDS-PAGE)
Sequence	MTKPITVQPG VQLRDADKMA LIPVKFLPEQ EGEQLKKPDW MRIRLPKTDE KIQNVKNIMR KNNLHSVCEE ASCPNLSECF NHGTATFMIL GAICTRHCPF CDVAHGKPLA PDADEPKKLA NTIREMALKY VVITSVDRDD LRDGGAQHFA DCIREIRLAS PNTRIETLTP DFRGRMDKAL DVFRETTPDV FNHNLETAPR LYSMARPGAD YAWSLKLQK MKELHPDLPT KSGLMMLGLGE TNDEIVQVLK DLRAHGVTML TLGQYLQPSR HHLPVKRYVP PQEFDELKAI ALDLGFTHAA CGPFVRRSSYH ADLQAQQQEV K
Source	Yeast
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.