



Recombinant Human Fibroblast growth factor 12 (FGF12)

Product Code	CSB-BP008618HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61328
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MESKEPQLKGIVTRLFSQQGYFLQMHPDGTIDGTKDENS DYTLFNLIPVGLRV VAIQGVKASLYVAMNGEGYLYSSDVFTPECKFKESV FENYYYVIYSSTLYRQQE SGRAWFLGLNKEGQIMKGNRVKKT KPSSHVFPKPIEVC MYREPSLHEIGE KQ GRSRKSSGTP TMNGGKVVNQDST
Source	Baculovirus
Target Names	FGF12
Protein Names	Recommended name: Fibroblast growth factor 12 Short name= FGF-12 Alternative name(s): Fibroblast growth factor homologous factor 1 Short name= FHF-1 Myocyte-activating factor
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length of Isoform 2
Target Details	This protein is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth, and invasion. This growth factor lacks the N-terminal signal sequence present in most of the FGF family members, but it contains clusters of basic residues that have been demonstrated to act as a nuclear localization signal. When transfected into mammalian cells, this protein accumulated in the nucleus, but was not secreted. The specific function of this gene has not yet been determined. Two alternatively spliced transcript variants encoding distinct isoforms have been reported.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.