



# Recombinant *Candida albicans* Methionine aminopeptidase 2 (MAP2)

<b>Product Code</b>	CSB-EP013716CZE-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	C4YSA9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Candida albicans</i> (strain WO-1) (Yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGATEGEDT KVIESKINEL NIDKPKLEDN NEAKGNGNGN ESGGDDDDDK EDDDDNDEIT EPSTSTSSGD KKKKKKNKNKK KKKKKIVSID SSYPEGIFPE GQWMEYPLED INSYRTTSEE KRYLDRQQNN KWQDFRKGAE IHRRVRHKAQ SSIRPGMTMI EIANLIEDSV RNYSGNDHTL KAGIGFPTGL SLNHVAAHYT PNTGDKLILK KDDIMKVDIG VHVNGRICDS AFTMTFNEDG KYDTIMQAVK EATYTGIKES GIDVRLNDIG AAIQEVMESE EMEENGKTYP IKCIKLNNGH NIDDFVIHSG KSVPIIANGD MTKMEEGETF AIETFGSTGN GYVLPEGECS HYAMNKGVEH LKPPSERSKQ LLETIKQNFG TLPWCRRYLE RTGEEKYLFA LNQLVRHGIV EEYPPIVDKR GSYTAQFEHT ILLHPHKKEV VTKGDDY
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	MAP2
<b>Protein Names</b>	Recommended name: Methionine aminopeptidase 2 Short name= MetAP 2 EC= 3.4.11.18 Alternative name(s): Peptidase M 2
<b>Expression Region</b>	1-447
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.