



Recombinant Human N-alpha-acetyltransferase 20 (NAA20)

| | |
|--------------------------|---|
| Product Code | CSB-EP015472HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P61599 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MTTLRAFTCD DLFRFNNINL DPLTETYGIP FYLQYLAHWP EYFIVAEAPG GELMGYIMGK AEGSVAREEW HGHVTSLSVA PEFRRRLGLAA KLMELLEIS ERKGGFFVDL FVRVSNQVAV NMYKQLGYSV YRTVIEYSSA SNGEPDEDAY DMRKALSRDT EKKSIIPLPH PVRPEDIE |
| Source | E.coli |
| Target Names | NAA20 |
| Protein Names | Recommended name: N-alpha-acetyltransferase 20 EC= 2.3.1.88 Alternative name(s): N-acetyltransferase 5 N-terminal acetyltransferase B complex catalytic subunit NAA20 N-terminal acetyltransferase B complex catalytic subunit NAT5 |
| Expression Region | 1-178 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |