



Recombinant Lithobates catesbeiana Ubiquitin-conjugating enzyme E2 S (ube2s)

Product Code	CSB-MP025477LQA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C1C3R6
Product Type	Recombinant Protein
Immunogen Species	Lithobates catesbeiana (American bullfrog) (Rana catesbeiana)
Purity	>85% (SDS-PAGE)
Sequence	MNSNVENLPP HVIRQVYKEV STLTDPPEG IKIIPNEEDI TDVQVSIIEGP EGTPYAGGIF RMKLILGKDF PAAPPKGYFL TKIFHPNVSS NGEICVNVLK KDWKAELGIR HVLLTIKCLL IHPNPESALN EEAGRLLLEN YEEYASRARL MTEIHAHSSS LRGKDPTDPC SSASVTGALG DGPMAKKHAG DRDKKLAACK KTDKKRALRR L
Source	Mammalian cell
Target Names	ube2s
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 S EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein S Ubiquitin-protein ligase S
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.