



Recombinant Brevibacillus brevis Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033BWL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C0Z4E2
Product Type	Recombinant Protein
Immunogen Species	Brevibacillus brevis (strain 47 / JCM 6285 / NBRC 100599)
Purity	>85% (SDS-PAGE)
Sequence	MGKLLPGLI HRLPQNAMSR TMGKITATPF SRLAIQRYIK HYQIDTSIIE KPASEYRTLK EFFSRRLKPA ARPIAPGPD IVSPVDGTVS QLGDICEGTL IQAKGKDFSV SELLGGSEEE AKRYYGKFI TIYLSRPDYH RIHMPVTGDL SSYCYLPGRL YPVNKLGIEN VDRLFARNER LVTHIKTDSL GDMALVKVGA LFVGSVKVCY NTATTNIKHG RQTHEKIAGT PRYEKGSELG WFEFG
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.