



Recombinant Mouse Abl interactor 2 (Abi2)

Product Code	CSB-BP001102MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P62484
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAELQMLLEE EIPGGRRALF DSYTNLERVA DYCENNYIQS PDKQRALEET KAYTTQSLAS VAYLINTLAN NVLQMLDIQA SQLRRMESSI NHISQTVDIH KEKVARREIG ILTTNKNTSR THKIIAPANL ERPVRYIRKP IDYTILDDIG HGKVVSTQNM KMGGLPRTP PTQKPPSPPM SGKGTLGRHS PYRTLEPVRP PVVPNDYVPS PTRNMAPSQQ SPVRTASVNQ RNRTYSSSGS SGGSHPSRS SSRENSGSGS VGVPIAVPTP SPPSVFPGHP VQFYSMNRPA SRHTPPTIGG SLPYRRPSI TSQTSLQNQM NGGPFYNQNP VSDTPPPPPP VEEPVFDESP PPPPPEDYE EEEAAVVEYS DPYAEEDPPW APRAYLEKVV AIYDYTKDKE DELSFQEGAI IYVIKKNDDG WYEGVMNGVT GLFPGNYVES IMHYSE
Source	Baculovirus
Target Names	Abi2
Protein Names	Recommended name: Abl interactor 2 Alternative name(s): Abelson interactor 2 Short name= Abi-2
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.