



# Recombinant Mouse Abl interactor 2 (Abi2)

<b>Product Code</b>	CSB-EP001102MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P62484
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAELQMLLEE EIPGGRRALF DSYTNLERVA DYCENNYIQS PDKQRALEET KAYTTQSLAS VAYLINTLAN NVLQMLDIQA SQLRRMESSI NHISQTVDIH KEKVARREIG ILTTNKNTSR THKIIAPANL ERPVRYIRKP IDYTIILDDIG HGKVVSTQNM KMGGLPRTP PTQKPPSPPM SGKGTLGRHS PYRTLEPVRP PVVPNDYVPS PTRNMAPSQQ SPVRTASVNQ RNRTYSSSGS SGGSHPSRS SSRENSGSGS VGVPIAVPTP SPPSVFPGHP VQFYSMNRPA SRHTPPTIGG SLPYRRPSI TSQTSLQNM NGGPFYNQNP VSDTPPPPPP VEEPVFDESP PPPPPEDYE EEEAAVVEYS DPYAEEDPPW APRAYLEKVV AIYDYTKDKE DELSFQEGAI IYVIKKNDDG WYEGVMNGVT GLFPGNYVES IMHYSE
<b>Source</b>	E.coli
<b>Target Names</b>	Abi2
<b>Protein Names</b>	Recommended name: Abl interactor 2 Alternative name(s): Abelson interactor 2 Short name= Abi-2
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.