



Recombinant *Micrococcus luteus* Lipoyl synthase (lipA)

Product Code	CSB-EP012927MTP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C5CCN0
Product Type	Recombinant Protein
Immunogen Species	<i>Micrococcus luteus</i> (strain ATCC 4698 / DSM 20030 / JCM 1464 / NBRC 3333 / NCIMB 9278 / NCTC 2665 / VKM Ac-2230) (<i>Micrococcus lysodeikticus</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSLAPEGRRLLRVEARNAEV PVERKPEWIK AKVHMGPEYI GLKNKVKASAG LHTVCEEAGC PNIFECWEDR EATFLIGGDI CTRRCDFCDI TSGKPRPLDM EEPQKVAENI REMDLRYATV TGVARDDQKD GAAWLYAETI RRIHALNPGT GVEILPPDFG AVPELVQQVF DARPEVFAHN LETVPRIFKR IRPAFTYEKS LRVLTMAKAD GLVTKSNLIL GMGEEDHEID QALVDLHESG CDIITITQYV RPSKLHHPID RWVKPQEFVQ WSQRAEEIGF QGVMAGPLVR SSYRAGKLYA QAMQRLGRTL PENLAHLAGE KTARQEASAV VAQMS
Source	E.coli
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.