



Recombinant *Drosophila ananassae* Ubiquitin-like modifier-activating enzyme 5 (GF23479)

Product Code	CSB-YP025415DKY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	D1GY43
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSHAIDELQA IIAELKVELE EQKTSNRQAR SRIDRMSAEV VDSNPYSRLM ALQRMNIVKE YERIRDKAVA VVGVGGVGSV TADMLTRCGI GKLILFDYDK VELANMNRLF FTPDQAGLSK VEAAARTLNF INPDVQIETH NYNITTVDNF DRFLATITES GKELGQPVDL VLSCVDNFEA RMAINAACNE RNLNWFESGV SENAVSGHIQ FIRPGDTACF ACAPPLVVAE NIDEKTLKRE GVCAASLPTT MGITAGFLVQ NALKYLLNFG EVSDYLG YNA LNDFPRMTL KPNPQCDDR N CLLRQKEFQA RPKPIEVKED VSSSDEPLHA TNEWGIELVA DDEPVNCPEP AKSSAVVQGL KLAYEAPKE KVEEENVATV SDETSLEDLM AQMKSM
Source	Yeast
Target Names	GF23479
Protein Names	Recommended name: Ubiquitin-like modifier-activating enzyme 5 Short name= Ubiquitin-activating enzyme 5
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.