



Recombinant *Saccharomyces cerevisiae* Methylthioribose-1-phosphate isomerase (MRI1)

Product Code	CSB-EP014799SVN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C8ZJE0
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SLEAIVFDR SEPENVSVKV LDQLLLPTYT KYVPIHTIDD GYSVIKSMQV RGAPAI AIVG SLSVLTEVQL IKHNPTSDVA TLYSLVNWES TKTVLNKRDL FLLSSRPTAV NLSNSLVEIK NILKSSSDLK AFDGSLYNYV CELIDEDLAN NMKMGDNGAK YLIDVLQKDG FKDEFVALTI CNTGSLATSG YGTALGVIRS LWKDSLAKTD KADSGLDNEK CPRMGHVFPL ETRPYNQGSR LTAYELVYDK IPSTLITDSS IAYRIRTSPI PIKAAFVGAD RIVRNGDTAN KIGTLQLAVI CKQFGIKFFV VAPKTTIDNV TETGDDIIVE ERNPEEFKVV TGTVINPENG SLILNESGEP ITGKVGIAPL EINVWNPAPD ITPHELIDGI ITEEGVFTKN SSGEFQLESL F
Source	E.coli
Target Names	MRI1
Protein Names	Recommended name: Methylthioribose-1-phosphate isomerase Short name= M1Pi Short name= MTR-1-P isomerase EC= 5.3.1.23 Alternative name(s): S-methyl-5-thioribose-1-phosphate isomerase Translation initiation factor eIF-2B sub
Expression Region	2-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.