



Recombinant Schizosaccharomyces pombe Exosome complex component rrp40 (rrp40)

Product Code	CSB-EP007891SXV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13903
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAEVEEKGLL SYYFPGERIP DSII SQD GSI RLG PGLV FQK KDREEEIVVS KAGRLHQ T GK NAMLIDSR TK RYIPATNDQV IGQIISRFAE GYRVDIGSAH IAQLNALAFE NVTRKSRPNL NVGSLVYARV SLADRDMEPE LECFDATTGK AAGYGELKNG YMITGLSLSH CRKLILPKNT LLQTLGSYIP FEIAVGMNGR VWVNSENLST TVLICTAIRN CEFMSDEEQI KYCKDLIKKL
Source	E.coli
Target Names	rrp40
Protein Names	Recommended name: Exosome complex component rrp40 Alternative name(s): Ribosomal RNA-processing protein 40
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.