



Recombinant Rat Cyclin-dependent kinase 18 (Cdk18)

Product Code	CSB-MP017650RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35832
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MNKMKNFKRR LLSVPRPET IEESLTFE TE QFNQLHTQRN EDGRDEPGQL SPGVQYQQRQ NQRFSMEDL NKRLSLPMDI RLPQEFLQKL QLENPGLPKP LTRMSRRASL SDIGFGKLET YVKLDKLGEG TYATVFKGRS KLTENLVALK EIRLEHEEGA PCTAIREVSL LKDLKHANIV TLHDLIHTDR SLTLVFEYLD SDLKQYLDHC GNLMNMHNVK IFMFQLLRGL AYCHRRKILH RDLKPQNLLI NERGELKLAD FGLARAKSVP TKTYSNEVVT LWYRPPDVLL GSTEYSTPID MWGVCILYE MATGKPLFPG STVKEELHLI FRLGTPTEE SWPGVTSISE FRAYNFPRYL PQPLLSHAPR LDTEGINLLT SLLLYESKSR MSAEAALSHP YFQSLGERVH QLDDTASIFS LKEIQLQKDP GYRGLAFQHP GRGKSRRQSI F
Source	Mammalian cell
Target Names	Cdk18
Protein Names	Recommended name: Cyclin-dependent kinase 18 EC= 2.7.11.22 Alternative name(s): Cell division protein kinase 18 PCTAIRE-motif protein kinase 3 Serine/threonine-protein kinase PCTAIRE-3
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.