



Recombinant Human Insulin-like growth factor 2 mRNA-binding protein 3 (IGF2BP3)

Product Code	CSB-MP011092HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00425
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNKLYIGNLS ENAAPSDLES IFKDAKIPVS GPFLVKTGYA FVDCPDESWA LKAIEALSGK IELHGKPIEV EHSVPKRQRI RKLQIRNIPP HLQWEVLDSL LVQYGVVESC EQVNTDSETA VVNVTYSSKD QARQALDKLN GFQLENFTLK VAYIPDEMAA QQNPLQQPRG RRGLGQRGSS RQGSPGSVSK QKPCDLPLRL LVPTQFVGAI IGKEGATIRN ITKQTQSKID VHRKENAGAA EKSITILSTP EG TSAACKSI LEIMHKEAQD IKFTEEIPLK ILAHNNFVGR LIGKEGRNLK KIEQD TDTKI TISPLQELTL YNPERTITVK GNVETCAKAE EEIMKKIRES YENDIASMNL QAHLIPGLNL NALGLFPPTS GMPPTSGPP SAMTPPYQF EQSETETVHL FIPALSVGAI IGKQGQHIKQ LSRFAGASIK IAPAEAPDAK VRMVIITGPP EAQFKAQGRI YGKIKEENFV SPKEEVKLEA HIRVPSFAAG RVIGKGGKTV NELQNLSSAE VVVPRDQTPD ENDQVVVKIT GHFYACQVAQ RKIQEILTQV KQHQQQKALQ SGPPQSRRK
Source	Mammalian cell
Target Names	IGF2BP3
Protein Names	Recommended name: Insulin-like growth factor 2 mRNA-binding protein 3 Short name= IGF2 mRNA-binding protein 3 Short name= IMP-3 Alternative name(s): IGF-II mRNA-binding protein 3 KH domain-containing protein overexpressed in canc
Expression Region	1-579
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.