



# Recombinant Human Nucleolar protein 56 (NOP56)

<b>Product Code</b>	CSB-BP015939HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O00567
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVLLHVLFEH AVGYALLALK EVEEISLLQP QVEESVLNLG KFHSIVRLVA FCPFASSQVA LENANAVSEG VVHEDLRLLL ETHLPSK KKK VLLGVGDPKI GAAIQEELGY NCQTGGVIAE ILRGVRLHFH NLVKGLTDLS ACKAQLGLGH SYSLAKVKFN VNRVDNMIQ SISLLDQLDK DINTFSMRVR EWYGYHFP EL VKIINDNATY CRLAQFIGNR RELNEDKLEK LEELTMDGAK AKAILDASRS SMGMDISAID LINIESFSSR VVSLSEYRQS LHTYLRKMS QVAPSLALI GEAVGARLIA HAGSLTNLAK YPASTVQILG AEKALFRALK TRGNTPKYGL IFHSTFIGRA AAKNKGRISR YLANKCSIAS RIDCFSEVPT SVFGEKLREQ VEERLSFYET GEIPRKNLDV MKEAMVQAE AAEITRKLE KQEKRLKKE KKRLAALALA SSENSSTPE ECEEMSEKPK KKKKQKPQEV PQENGMEDPS ISFSKPKKKK SFSKEELMSS DLEETAGSTS IPKRKKSTPK EETVNDPEEA GHRSGSKKKR KFSKEEPVSS GPAAVAVGKSS SKKKKKFHKA SQED
<b>Source</b>	Baculovirus
<b>Target Names</b>	NOP56
<b>Protein Names</b>	Recommended name: Nucleolar protein 56 Alternative name(s): Nucleolar protein 5A
<b>Expression Region</b>	1-594
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.