



# Recombinant Drosophila melanogaster Eukaryotic translation initiation factor 3 subunit I (Trip1)

<b>Product Code</b>	CSB-BP007538DLU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O02195
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRPLMLQGHE RSITQIKYNR EGDLLFSCSK DQKPNVWYSL NGERLGTYDG HQGAVWCLDV DWESRKLITG AGDMTAKIWD VEYGTVIASI PTKSSVRTSN FSFSGNQAAY STDKAMGQSC ELFLIDVRNA DSSLSEQEPT LRIPMTESKI TSMLWGPLDE TIITGHDNGN IAIWDIRKGG KVVDSGTDHS AGINDMQLSK DGTMFVTASK DTTAKLFDSE SLMCLKTYKT ERPVNSAAIS PIMDHVVLGG GQDAMEVTTT STKAGKFDSR FFHLYEEEF ARLKGHFGPI NSLAFHPDGK SYASGGEDGF VRVQTFDSTY FENIFE
<b>Source</b>	Baculovirus
<b>Target Names</b>	eIF3i
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit I Short name= eIF3i Alternative name(s): Eukaryotic translation initiation factor 3 subunit 2 TRIP-1 homolog eIF-3-beta
<b>Expression Region</b>	1-326
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.