



# Recombinant *Drosophila melanogaster* UDP-glucose 6-dehydrogenase (sgl)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP025567DLU-B  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | O02373   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | <i>Drosophila melanogaster</i> (Fruit fly)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MKVCCIGAGY VGGPTCAVMA LKCPDIVITL VDKSSERIAQ WNSDKLPIYE<br>PGLDEVVKKC RNVNLFSTD IETAIKEADL IFISVNTPTK TCGNGKGRAA<br>DLKYVESAAR MIAEIAQSNK IVVEKSTVPV RAAESIMHIL RANQKPGIHY<br>DILSNPEFLA EGTAINDLLN ADRVIGGEE TPEGHQAWEK LSWIYEHWIP<br>KQNILTTNTW SSELSKLAAN AFLAQRISI NSLSAVCEAT GADVSEVARA<br>VGLDSRIGSK FLQASVGFEGG SCFQKDILNL IYICENLNLP EVAAYWQQVI<br>DMNEYQKRRF SQKIIESLFN TVSDKRIAIL GFAPKNTGD TRETAAITVC<br>QTLLEEGAAL DIYDPKVEPE QIDDLTHPS VTESPEKVKK AVQIHSDPYS<br>AVRATHALVI CTEWDEFVDL DFKRIYQSMK KPAYIFDGRK ILDBERLQKI<br>GFHVQTIGKK YQRTGLLRSW GIVPQL |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | sgl  |
| <b>Protein Names</b>     | Recommended name: UDP-glucose 6-dehydrogenase Short name= UDP-Glc dehydrogenase Short name= UDP-GlcDH Short name= UDPGDH EC= 1.1.1.22<br>Alternative name(s): Protein sugarless Protein suppenkasper   |
| <b>Expression Region</b> | 1-476  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |