



Recombinant Arabidopsis thaliana Eukaryotic translation initiation factor 3 subunit F (TIF3F1)

Product Code	CSB-MP007535DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O04202
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAATSEHTIL QFVSPSSTAS ATTSVLTARI HPLVIFNVCD CFVRRPDSAE RVIGTLLGSI LPDGTVDIRN SYAVPHNESS DQVAVDIDYH HNMLASHLKV NSKETIVGWY STGAGVNGGS SLIHDFYARE VPNIHLTVD TGFTNGEGTI KAFVSSNLSL GDRQLVAHFQ EIPVDLRMVD AERVGFDVLK ATSVDKLPND LEGMELTMER LLTLINDVYK YVDSVVGQI APDNNIGRFI ADAVASLPKL PPQVFDNLVN DSLQDQLLLL YLSSITRTQL SLAEKLNTAA QML
Source	Mammalian cell
Target Names	TIF3F1
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit F Short name= eIF3f Alternative name(s): eIF-3-epsilon eIF3 p32 subunit
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.