



Recombinant Schizosaccharomyces pombe V-type proton ATPase subunit H (vma13)

Product Code	CSB-MP002413SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14265
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSNSDLELSN AASPPPVELD NSQVDEIINN VRCVAIPWQG YQRSGSLEEN ELQEIENLTG KPLSAYVKTA EEDTTAYSNL FLKLLSMKDT PDVNFALVK LADTLLNSNK FLSAFGPAFY DFLEKDESYI NYLDDDSKLL FARVFALCSS SSPCSVAKAF TLFLEYLGKL MQSLNPLTRL FAVQCLNGVL TLKAHRYALW AENTCSFRLA ELLRNSIGDT QLQYYSLFCF WQLTFESHIA QDINKRFDLI KLLVQIIRSD TKTKVYRLVL AILVNLIDKA PKDTISTMLL EHVDKAVQLL QKRKWADEDI TNYLDFITST LDESSKHLST FDMYKSELDT GILHWSPSHR SEDFWHQNAK RLNEDNYALL KKLHFHIVQYN EDNTSLAVAC HDLGAYIRSY PEGRSLIIKY GAKQRIMDLM SHPDPEVRFE ALSTVQLLMT EVCFLSKITL
Source	Mammalian cell
Target Names	vma13
Protein Names	Recommended name: V-type proton ATPase subunit H Short name= V-ATPase subunit H Alternative name(s): V-ATPase 54 kDa subunit Vacuolar proton pump subunit H
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.