



MDK Antibody, Biotin conjugated

Product Code	CSB-PA10574D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21741
Immunogen	Recombinant Human Midkine protein (21-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Midkine (MK) (Amphiregulin-associated protein) (ARAP) (Midgestation and kidney protein) (Neurite outgrowth-promoting factor 2) (Neurite outgrowth-promoting protein), MDK, MK1 NEGF2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MDK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.