



pstS1 Antibody, Biotin conjugated

Product Code	CSB-PA14919D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WGU1
Immunogen	Recombinant Mycobacterium tuberculosis Phosphate-binding protein PstS 1 protein (24-373AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphate-binding protein PstS 1 (PBP 1) (PhoS1) (PstS-1) (38-kDa glycolipoprotein) (38-kDa lipoprotein) (P38) (Antigen Ag78) (Protein antigen B) (Pab), pstS1, phoS1
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Research Area	Others
Target Names	pstS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.