




Hpx Antibody

Product Code	CSB-PA010723LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20059
Immunogen	Recombinant Rat Hemopexin protein (24-460AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Hpx

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate(30µg), Mouse Brain tissue lysate(30µg), Mouse Liver tissue lysate(30µg), Mouse Spleen tissue lysate(30µg), Mouse Kidney tissue lysate(30µg), Rat Brain tissue lysate(30µg), Rat Liver tissue lysate(30µg), Rat Kidney tissue lysate(30µg)

All lanes: Hpx antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/40000 dilution

Predicted band size: 52.0 kDa

Observed band size: 52.0 kDa

Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.