



SPA17 Antibody, FITC conjugated

Product Code	CSB-PA12124C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15506
Immunogen	Recombinant Human Sperm surface protein Sp17 protein (1-146AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm surface protein Sp17 (Cancer/testis antigen 22) (CT22) (Sp17-1) (Sperm autoantigenic protein 17) (Sperm protein 17), SPA17, SP17
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPA17