



SPA17 Antibody, FITC conjugated

Product Code	CSB-PA12124C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15506
Immunogen	Recombinant Human Sperm surface protein Sp17 protein (1-146AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm surface protein Sp17 (Cancer/testis antigen 22) (CT22) (Sp17-1) (Sperm autoantigenic protein 17) (Sperm protein 17), SPA17, SP17
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPA17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.