



CASP3 Antibody, FITC conjugated

Product Code	CSB-PA05689C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42574
Immunogen	Recombinant Human Caspase-3 protein (29-175AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caspase-3 (CASP-3) (EC 3.4.22.56) (Apopain) (Cysteine protease CPP32) (CPP-32) (Protein Yama) (SREBP cleavage activity 1) (SCA-1) [Cleaved into: Caspase-3 subunit p17; Caspase-3 subunit p12], CASP3, CPP32
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CASP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.