



Corticosterone Monoclonal Antibody

Product Code	CSB-MA004311I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Corticosterone-BSA conjugate
Raised In	mouse
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	11-dehydro-17-hydroxycorticosterone, cortisone
Product Type	Small Molecules Antibodies
Target Names	CORT
Accession NO.	4G1B4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.