



Recombinant Protein LTV1 homolog (T23D8.3)

Product Code	CSB-BP013231CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02327
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGKRKAFIDK KASQSFRLVP DNRERSERFK PTQEHLEEQQ KYGVYYDDDY DYLQHMRAVN EPMKLENVHE EVEKTTIKTS GSSAFPPAPP LFGLVGALKK PEFFDEDVAN ALEEVTDDRN TGELEDNFIT LAGLLDERT TVYRSTRRGE DSEEEEDDDE DDEMYDDYND DELFGEEAVG EIRVERADQR VIDNAFEELM DREYNTDQIG ELDGDDYDVG GALEPNAGRL HKLINDKGPS NAEYDEELAK HYVRERMRLI EEGVIKDKEE YEIVEVDEGT NKKMKWDCES FATQYTNIYN HPTLIKEPRG LSRKALKRFD KAVEEMDIAE EDEDDDEDME DDDDKESVFS TVSTFRPKNE TPEQRSLRKK AVKEARKLRR VEKKANKTMF AEEKRKLAKG RIGQIKARPI
Source	Baculovirus
Target Names	T23D8.3
Protein Names	Recommended name: Protein LTV1 homolog
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.