



CDKN1A Antibody, FITC conjugated

Product Code	CSB-PA005086DC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38936
Immunogen	Recombinant Human Cyclin-dependent kinase inhibitor 1 protein (2-164AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase inhibitor 1 (CDK-interacting protein 1) (Melanoma differentiation-associated protein 6) (MDA-6) (p21), CDKN1A, CAP20 CDKN1 CIP1 MDA6 PIC1 SDI1 WAF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDKN1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.