



MSA2 Antibody

Product Code	CSB-PA343481DA01EWP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50498
Immunogen	Recombinant Plasmodium falciparum Merozoite surface antigen 2 protein (109-246AA)
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Merozoite surface antigen 2 (MSA-2) (45 kDa merozoite surface antigen), MSA2
Immunogen Species	Plasmodium falciparum
Research Area	Others
Target Names	MSA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.