



KLK3 Antibody, Biotin conjugated

Product Code	CSB-PA06779D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07288
Immunogen	Recombinant Human Prostate-specific antigen protein (28-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostate-specific antigen (PSA) (EC 3.4.21.77) (Gamma-seminoprotein) (Seminin) (Kallikrein-3) (P-30 antigen) (Semenogelase), KLK3, APS
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KLK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.