



ANXA6 Antibody, HRP conjugated

Product Code	CSB-PA00794B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08133
Immunogen	Recombinant Human Annexin A6 protein (2-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A6 (67 kDa calelectrin) (Annexin VI) (Annexin-6) (Calphobindin-II) (CPB-II) (Chromobindin-20) (Lipocortin VI) (Protein III) (p68) (p70), ANXA6, ANX6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ANXA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.