



HDGF Antibody, Biotin conjugated

Product Code	CSB-PA01004D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51858
Immunogen	Recombinant Human Hepatoma-derived growth factor protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatoma-derived growth factor (HDGF) (High mobility group protein 1-like 2) (HMG-1L2), HDGF, HMG1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	HDGF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.