



# Recombinant Human BPI fold-containing family B member 3 (BPIFB3)

<b>Product Code</b>	CSB-BP003693HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P59826
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	CQELLETVGT LARIDKDELG KAIQNSLVGE PILQNVLGSV TAVNRGLLGS GGLLGGGGLL GHGGVFGVVE ELSGLKIEEL TLPKVLLKLL PGFGVQLSLH TKVGMHCSGP LGLLQLAAE VNVTSRVALA VSSRGTPILI LKRCSTLLGH ISLFSGLLPT PLFGVVEQML FKVLPGLLCP VVDSVLGVVN ELLGAVLGLV SLGALGSVEF SLATLPLISN QYIELDINPI VKSVAGDIID FPKSRAPAKV PPKGDHTSQV MVPLYLFNTT FGLLQTNGAL DMDITPELVP SDVPLTTTDL AALLPEALGK LPLHQQLLLF LRVREAPT VT LHNKKALVSL PANIHVLFYV PKGTPESLFE LNSVMTVRAQ LAPSATKLHI SLSLERLSVK VASSFTHAFD GSRLEEWLSH VVGAVYAPKL NVALDVGIPL PKVLNINFSN SVLEIVENAV VLTVAS
<b>Source</b>	Baculovirus
<b>Target Names</b>	BPIFB3
<b>Protein Names</b>	Recommended name: BPI fold-containing family B member 3 Alternative name(s): Ligand-binding protein RYA3 Long palate, lung and nasal epithelium carcinoma-associated protein 3
<b>Expression Region</b>	21-476
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.