



Recombinant Human BPI fold-containing family B member 3 (BPIFB3)

Product Code	CSB-EP003693HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59826
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	CQELLETVGT LARIDKDELG KAIQNSLVGE PILQNVLGSV TAVNRGLLGS GGLLGGGGLL GHGGVFGVVE ELSGLKIEEL TLPKVLLKLL PGFGVQLSLH TKVGMHCSPG LGLLQLAAE VNVTSRVALA VSSRGTPIIL LKRCSTLLGH ISLFSGLLPT PLFGVVEQML FKVLPGLLCP VVDSVLGVVN ELLGAVLGLV SLGALGSVEF SLATLPLISN QYIELDINPI VKSVAGDIID FPKSRAPAKV PPKKDHTSQV MVPLYLFNTT FGLLQTNGAL DMDITPELVP SDVPLTTTDL AALLPEALGK LPLHQQLLLF LRVREAPTVT LHNKKALVSL PANIHVLFYV PKGTPESLFE LNSVMTVRAQ LAPSATKLHI SLSLERLSVK VASSFTHAFD GSRLEEWLSH VVGAVYAPKL NVALDVGIPL PKVLNINFSN SVLEIVENAV VLTVAS
Source	E.coli
Target Names	BPIFB3
Protein Names	Recommended name: BPI fold-containing family B member 3 Alternative name(s): Ligand-binding protein RYA3 Long palate, lung and nasal epithelium carcinoma-associated protein 3
Expression Region	21-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.