



ISG15 Antibody, HRP conjugated

Product Code	CSB-PA09744B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05161
Immunogen	Recombinant Human Ubiquitin-like protein ISG15 protein (2-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like protein ISG15 (Interferon-induced 15 kDa protein) (Interferon-induced 17 kDa protein) (IP17) (Ubiquitin cross-reactive protein) (hUCRP), ISG15, G1P2 UCRP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ISG15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.